

TC4-SOP 0-2b

ELISA titer test

Project Name: _____

MATERIALS:

- **Blocking Buffer**
Gelatin-NET
10% FBS in 1X PBS
- **Coating Buffer - 0.1 M Sodim Carbonate, pH 9.5**
7.13g NaHCO₃, 1.59g Na₂CO₃; q.s. to 1.0 L; pH to 9.5 with 10 N NaOH. Freshly prepare or use within 7 days of preparation, stored at 2-8°C
- **PBST buffer**
PBS buffer with 0.05% (v/v) Tween-20
- **Test antibody samples (1st antibody)**
- **2nd antibody (HRP-goat-anti-mouse)**
- **TMB substrate - Tetramethylbenzidine (TMB) and Hydrogen Peroxide**
KPL SureBlue™ TMB Microwell Peroxidase Substrate (Cat. No. 52-00-00).
BD Phamingen™ TMB Substrate Reagent Set (Cat. No. 555214).
- **Stop Solution**
TMB – 1 M H₃PO₄ or 2 N H₂SO₄.
- **ELISA plate (NUNC 442404)**
- **12-channel pipette**
- **ELISA reader**

METHODS:

- Coat plate with antigen :**
 1. Prepare an antigen solution in Coating Buffer at 0.2~10 µg/ml. The concentration of antigen is usually 10 µg/ml. Prepare ~10 ml antigen solution for each plate.
 2. Using a 12-channel pipette and tips, dispense 100 µL antigen solution into each well of ELISA plate.
 3. Seal coated plates using 96-well adhesive plastic seal and incubate overnight at room temperature (RT) or 2 hr at 37°C .
 - Prepare control: **NC**-Coating buffer without antigen ; **PC**-Coating buffer with sera.
- Block :** Aspirate wells and wash 3 times PBST buffer. After last wash, invert plate and blot on absorbent paper to remove any residual buffer. Block plates with ≥ 200 µL/well blocking buffer. Incubate at RT for 1 hour. Aspirate and wash 3 times with PBST.
- Add 1st Antibody:** Pipette 100 µL of each sample and control diluted in blocking buffer into appropriate coated wells. Seal plate and incubate ≥ 1 hr at RT. Aspirate and wash 3 times with PBST.
- Add 2nd Antibody :** Pipette 100 µL of 5000X 2nd Antibody diluted in blocking buffer to each well. Seal plate and incubate for 1 hr at RT. Aspirate and wash 7 times with PBST buffer.
- Add TMB substrate :** Add 100 µL of TMB substrate solution to each well. Incubate plate (without plate sealer) for 10~30 min at RT in the dark and then add 50 µL of Stop solution to each well.
- Optical density :** Using ELISA reader to read absorbance at 450 nm within 30 min of stopping reaction. If wavelength correction is available, subtract absorbance 570 nm from absorbance 450 nm.

Note:

Date	Operator	QC