

TC4-SOP 1-1m

 Immunize mice

Project Name: _____

MATERIALS:

- Mouse x _____
 - BALB/c Birth Date: ____/____/____
- 3-way stopcock (Nipro)
- Antigen: _____
 - _____ $\mu\text{g}/\mu\text{L}$, _____ μL _____ μg
 - Powder Gel liquid, solvent _____
- Adjuvant
 - Complete Freund's Adjuvant (CFA) Incomplete Freund's Adjuvant (IFA)
 - other: _____ None
- Sterile syringes (Terumo, 27 GX1/2", 1 mL)
- Paraffin membrane (Parafilm)
- PBS buffer (1X)

METHODS:

- For first immunization, prepare antigen using 100 μg protein or peptide per animal to be immunized. For subsequent boosters, using 50 μg antigen. Add 1X PBS to 250 μL . (Antigen _____ μL , add 1X PBS to 250 μL)
- Add antigen into a sterile syringe using tips. Carefully clean the end of syringe without needle, connect syringe to a 3-way stopcock. Seal the connection part by using parafilm.
- Add a volume of CFA¹/IFA² equal to the antigen volume into another syringe using tips. Carefully clean the end of syringe, connect syringe to the 3-way stopcock. Seal the connection part by using parafilm.
 1. Completely resuspend CFA to disperse the *Mycobacterium tuberculosis* bacilli.
 2. CFA is for first booster, IFA is for subsequent booster.
- Turn the 3-way stopcock, **discharging antigen into adjuvant***.
 - * Add one drop of adjuvant into the syringe plunger.
- Discharging back and forth half an hour until a thickened mixture results*.
 - * A stable emulsion will not disperse when a drop of it is placed in water, or looks like the tooth paster.
- Transfer all of the adjuvant/antigen emulsion to one syringe. Attach a sterile 27G needle to the syring containing the emulsion.
- Inject emulsion intraperitoneally into the animal using <500 $\mu\text{L}/\text{mouse}$.
- Boost animal after 14 days and bleed to monitor the antibody titer.**
 - * If cell fusion is planned for 3~7 days after boosting, immunize with antigen alone in 1X PBS. If a fusion is not immediately planned, boost the animal with antigen emulsified in IFA.

Reference:

John E. C., Barbara E. B., David H. M., Ethan M.S., and Warren S. (eds.) 2009. *Current Protocol in immunology*. John Wiley and Sons. Inc.

Date	Operator	QC